



# Bioanalyzer Applications for Next-Gen Sequencing: Updates and Tips

*March 1<sup>st</sup>, 2011*

***Charmian Cher, Ph.D***  
***Field Applications Scientist***





**Agilent Technologies**

# Agenda


- 1** Next-gen sequencing library preparation workflow and important quality control steps.
- 2** Introduction to the 2100 Bioanalyzer.
- 3** The use of Bioanalyzer assays for:
  - assessing quality of starting material
  - monitoring size distribution after fragmentation and adapter ligation
  - quantifying yield and detection of artifacts post-PCR amplification
  - detecting small quantities of DNA in amplification-free protocols
- 4** Best practice tips to always obtaining accurate and reproducible data when using the Bioanalyzer.

# Next-Gen Sequencing



## Roche/454 and Life Technologies



GS FLX Titanium



## Illumina



SOLiD 4



HiSeq 2000

Adapted from Metzker, M. (2009) Nature Reviews


# Quality of DNA libraries is critical for sequencing success

For every application of next-gen sequencing eg. genome sequencing, transcriptome sequencing (RNA-seq), chromatin immunoprecipitation sequencing (ChIP-seq) or targeted resequencing, there is a specific protocol to convert the source nucleic acid to standard DNA libraries.

The aim of recent developments in library preparation methods is to produce a high-quality representative, non-biased DNA library from small amounts of starting material.

# Quality Control during Library Preparation


## Genome sequencing




*Adapted from Quail et al (2008) Nature Methods*

# Quality Control during Library Preparation

## SureSelect™ Target Enrichment System




# Implementing quality control with the 2100 Bioanalyzer



**SureSelect Human All Exon Kit**  
Illumina Paired-End Sequencing Library Prep


**Protocol**  
Version 18, October 2009  
SureSelect platform manufactured with Agilent SurePrint Technology  
Research Use Only. Not for use in Diagnostic Procedures.



**Step 3. Assess quality with Agilent 2100 Bioanalyzer**


Use a Bioanalyzer DNA 1000 chip and reagent kit.

- 1 Check that the 2100 Bioanalyzer electrodes have been cleaned as instructed in the reagent kit guide.
- 2 Open the Agilent 2100 expert software (version B.02.02 or higher), turn on the 2100 Bioanalyzer and check communication.
- 3 Prepare the chip, samples and ladder as instructed in the reagent kit guide.
- 4 Load the prepared chip into the 2100 Bioanalyzer and start the run within five minutes after preparation.
- 5 Within the instrument context, choose the appropriate assay from the drop down list.
- 6 Start the run. Enter sample names and comments in the Data and Assay context.
- 7 Verify the results.  
Check that the electropherogram shows a distribution with a peak height of 150 ± 10% nucleotides.




**Figure 2** Analysis of sheared DNA using a DNA 1000 Bioanalyzer assay. The electropherogram shows a distribution with a peak size of 150 ± 10% nucleotides.

All major next-gen sequencing providers recommend using the Bioanalyzer at various steps in the library preparation workflow.




**Genome Sequencer System**  
Application Note No. 5 | February 2007

**Amplicon Sequencing**



www.roche-applied-science.com

**Roche**




Agilent's gel-free sequencing process can improve library preparation at a large number of steps, starting with adapter clean-up and the PCR. The new SureSelect protocol provides an improved PCR cleanup. It uses a novel bead-based purification process and a different chemistry, combination of PCR conditions and a cycle time of 100 seconds per reaction. These improvements result in a more consistent PCR cleanup and improved library quality.

**Preparing Samples for Analysis of Small RNA**  
FOR RESEARCH ONLY

Topics

- 3 Introduction
- 4 Kit Contents and Equipment Checklist
- 5 Isolate Small RNA by Denaturing PAGE
- 6 Ligate 5' RNA Adapters
- 7 Ligate 3' RNA Adapters
- 8 Reverse Transcribe and Amplify the Small RNA Ligated with Adapters
- 9 Purify the Amplified cDNA Construct
- 10 Validate the Library

© 2008 Agilent  
Part # 526227 Rev. A  
March 2008




**SOLID™ Small RNA Expression Kit**  
(Part Number 430762)  
Protocol

- I. Introduction ..... 5
  - A. Background and Product Description
  - B. Procedure Overview
  - C. Reagents Provided with the Kit and Storage Conditions
  - D. Materials Not Provided with the Kit
  - E. Related Products Available from Applied Biosystems
- II. Procedure ..... 7
  - A. RNA Sample Type and Name
  - B. Hybridization Ligases
  - C. Reverse Transcription and Ethanol Digestion
  - D. Small RNA Library Amplification
  - E. Amplified Small RNA Library Cleanup
  - F. Size Selection of Amplified Small RNA Library by PAGE
- III. Troubleshooting ..... 18
  - A. Positive Control Reaction
  - B. No PCR Products in the Expected Size Range
- IV. Appendix ..... 20
  - A. Sequence of the SOLiD™ PCR Primers Included in the Kit
  - B. Filter/Strand Gel Instructions
  - C. References
  - D. Quality Control
  - E. Safety Information

**Fragment the Genomic DNA**

This protocol fragments the genomic DNA using a nebulization technique, which fragments DNA to less than 800 bp in minutes using a cost-effective, disposable device. Nebulization generates double-stranded DNA fragments that are blunt-ended or comprised of 3' or 5' overhangs.



**Figure 3** Fragment Genomic DNA

**Validate the Library**

The amount of starting material is very low (10 ng), and after 18 cycles of PCR, you will have a low number of molecules. It is recommended that you use a sensitive qPCR measurement assay such as the Quant-iT 2.0 qPCR Assay Kit, 150 assays (0.2-100 ng for use with the Quant-iT fluorometer (Invitrogen)). Note that this will not allow you to check the size and purity of your sample. Do not use an OD260/280 ratio for concentration measurements, since this will not distinguish dDNA from primers, and therefore cannot be used to validate the library.

**Bioanalyzer Method**

1. Load 1 µl of the resuspended construct and 1 µl of the negative control on an Agilent Technologies 2100 Bioanalyzer.
2. Check the size, purity, and concentration of the sample.

**Alternative Methods**


For users who do not have access to an Agilent Technologies 2100 Bioanalyzer or similar instrument, you may try using a sensitive qDNA measurement assay such as the Quant-iT 2.0 qPCR Assay Kit, 150 assays (0.2-100 ng for use with the Quant-iT fluorometer (Invitrogen)). Note that this will not allow you to check the size and purity of your sample. Do not use an OD260/280 ratio for concentration measurements, since this will not distinguish dDNA from primers, and therefore cannot be used to validate the library.

# Agenda

- 1** Next-gen sequencing library preparation workflow and important quality control steps.
- 2** Introduction to the 2100 Bioanalyzer.
- 3** The use of Bioanalyzer assays for:
  - assessing quality of starting material
  - monitoring size distribution after fragmentation and adapter ligation
  - quantifying yield and detection of artifacts post-PCR amplification
  - detecting small quantities of DNA in amplification-free protocols
- 4** Best practice tips to always obtaining accurate and reproducible data when using the Bioanalyzer.



# 2100 Bioanalyzer



**Sample volumes 1 - 5  $\mu$ l**

**10 -12 samples depending on Assay**

**Separation, staining, detection of samples**

**Results in 5-30 minutes available**

**No extra waste removal needed**

**Disposable Chip, no cross-contamination**




# Principle of Electrodriven Flow

The sample moves from the sample well through the micro-channels


The sample is injected into the separation channel

Sample components are electrophoretically separated

Components are detected by their fluorescence and translated into gel-like images (bands) and electropherograms (peaks)



# Bioanalyzer Kit Portfolio



DNA1000

DNA7500

DNA12000

High Sensitivity DNA

## DNA Assays:

- Sizing
- Quantitation
- PCR products, digests, larger DNA fragments
- 12 samples in 30 min.



6000 Nano

6000 Pico

Small RNA

## RNA Assays:

- Quantitation (Sizing in Small RNA)
- total RNA, mRNA
- purity & integrity determination
- 10 samples in 30 min.



Flow Cytometry

## Cell Assays:

- Analysis of 6 samples
- Two color detection
- Analysis of protein expression in cells



P230

P80

HSP 250

## Protein Assays:

- Sizing
- Quantitation
- cell lysates, column fractions, purified proteins, antibodies etc.
- 10 samples in 40 min.




# Agenda


- 1** Next-gen sequencing library preparation workflow and important quality control steps.
- 2** Introduction to the 2100 Bioanalyzer.
- 3** The use of Bioanalyzer assays for:
  - assessing quality of starting material
  - monitoring size distribution after fragmentation and adapter ligation
  - quantifying yield and detection of artifacts post-PCR amplification
  - detecting small quantities of DNA in amplification-free protocols
- 4** Best practice tips to always obtaining accurate and reproducible data when using the Bioanalyzer.

# 1 Integrity of starting material

- When performing RNA-seq, it is necessary to assess the integrity of total RNA prior to mRNA or small RNA isolation using the RNA 6000 Nano Assay.
- A RNA Integrity Number (RIN) of at least 8 is the recommended threshold.
- The quality of chromatin immunoprecipitated (ChIP) DNA can be assessed using the High Sensitivity Assay.




# Total RNA Integrity



# Analyzing small amounts of CHIP DNA


**High Sensitivity DNA Assay : concentration range of 100 – 10000 pg/μl**



**ChIP Post-IP**

**Average Size: 827 bp**

**Conc: 335 pg/μl**





**ChIP Post-IP**

**High molecular weight DNA  
contamination**

- **Fragment DNA prior to IP using Covaris, sonication or enzymatic digestion.**

## 2 Monitoring size distribution after fragmentation



- Comparison of fragmentation by nebulization with AFA technology (Covaris).
- For a 200-bp ( $\pm 20$  bp) library, the yield produced by AFA was four- to fivefold greater than that produced by nebulization.

A large genome center's improvements to the Illumina sequencing system

Michael A Quail, Iwanka Kozarewa, Frances Smith, Aiswya Scally, Philip J Stephens, Richard Durbin, Harold Swerdlow & Daniel J Turner

NATURE METHODS | VOL. 5 NO. 12 | DECEMBER 2008 | 1005




# The importance of monitoring size distribution


- **Choice of fragmentation can significantly affect the recovery of desired fragments and hence the amount of starting material required.**
- **For example, the Covaris typically produces narrower fragment distributions than nebulization, resulting in 4-fold greater recovery.**
- **Sequencing fragments that do not fall within the recommended size distribution may lead to low read depth or even a lack of read coverage for specific portions of the sequence.**
- **When performing SureSelect target enrichment, fragment size distribution can affect final % on-target capture.**
- **Probe-based sonication methods routinely introduce sample-to-sample variability and are not recommended.**

# RNA Fragmentation

## SOLiD




Poly(A) RNA fragmented by RNase III




rRNA-depleted RNA fragmented by RNase III

## Roche/454




Non-fragmented mRNA



Fragmented mRNA


# DNA Fragmentation

DNA sheared using the Covaris in the SureSelect<sup>XT</sup> Illumina Paired-End sequencing protocol.



Size distribution with peak height between 150-200nt.



# DNA Fragmentation - Tailing



**Cause :** Too much DNA was used for shearing on Covaris, resulting in incomplete shearing. For example, Covaris settings in the SureSelect protocol are optimized to yield a peak size range of 150-200bp when starting with a maximum of 3 $\mu$ g gDNA. Adhere to recommended amounts of starting material.

*Maria-Celeste Ramirez, PhD  
Agilent Technologies*


# DNA Fragmentation – Uneven shearing



## Causes :

- Poor DNA quality.
- DNA eluted in wrong buffer (not TE).
- Covaris issues
  - Sample volume less than 120 $\mu$ l. This allows an air gap, resulting in inconsistent fragmentation.
  - Water level is too high/low.
  - Temperature of waterbath not between 6-8 $^{\circ}$ C.
  - Insufficient degassing.


# DNA Fragmentation – Increased size range



**Cause** : Covaris issues. Perform a control experiment using DNA of known quality, such as commercially available lambda DNA.


### 3 Monitoring size distribution after adapter-ligation

Post - fragmentation




Peak size of 190bp

After adaptor ligation




Peak size of 250 ± 10bp



- Adaptor ligation is not 100% efficient so the DNA 1000 assay can be used to visualize the size shift after successful ligation of adaptors.
- Inefficient adapter ligation will result in reduced library complexity after PCR.

# Adapter ligation – Excess adapters



**Cause** : Inefficient ligation due to too much input DNA or the use of incorrect ligation temperature (ligation is performed at 20-25°C. When using a PCR machine, make sure the lid is not heated).



# Automated size selection after adapter ligation

In whole genome sequencing, small-RNA seq and ChIP-seq protocols, a gel-based size selection step is performed after adapter ligation (SureSelect™ uses SPRI beads).

This can be automated using the Pippin Prep™ DNA size selection system, which uses pre-cast agarose gel cassettes and elutes size-selected fragments in buffer.



Size selection  
using Pippin Prep™




Quality control of fractions using  
the Bioanalyzer




Sequencing

# Automated size selection after adapter ligation

- Resolving Illumina ligation reactions directly on gel have been observed to be affected by the DNA-bound ligase, resulting in extraction of incorrectly-sized DNA.
- A simple protocol was developed that includes an additional proteinase K digestion for 10min at 37°C prior to gel electrophoresis.
- Samples were then loaded on a 2% Pippin Prep cassette using settings to collect DNA fragments of 335bp (range: 308-362bp). Extracted samples analyzed on Bioanalyzer.



No treatment



Proteinase K treatment


## Pippin Prep™ DNA Size Selection System

A Rapid Proteinase K Cleanup Method for Illumina Adapter Ligation Reactions Prior to Pippin Prep Fractionation

Matthew Mayho and Michael Quail. The Wellcome Trust Sanger Institute, Wellcome Trust Genome Campus, Hinxton, Cambridge CB10 1SA, UK.


## 4 Quality control after PCR amplification

- PCR amplification is necessary in library preparation to enrich for adapter-ligated fragments, as well as to add indexes.
- In the SureSelect protocol, PCR is also used to amplify captured DNA for final QC and quantification.
- PCR can create bias since some regions of template DNA amplify more poorly than others (eg, GC-rich regions). PCR artifacts caused by overamplification or primer dimers can also affect sequencing coverage and accuracy.




# PCR amplification – detecting artifacts

## High Quality DNA library



## Primer dimers and PCR artifact



Primer dimers

PCR artifact



DNA library (50% of total DNA)

# Dealing with primer dimers and artifacts


- **Primer dimers present in the library will be sequenced, resulting in wastage of sequencing capacity and increased duplicate reads.**
- **When working with small amounts of DNA, such as in ChIP-seq or SureSelect™ target enrichment, primer dimers can be removed by performing additional SPRI bead clean-up steps.**
- **Gel-based size selection can also be performed to remove primer dimers.**
- **Repeat PCR with fewer cycles to prevent formation of artifacts using remaining library.**

# Reducing number of amplification cycles

Illumina GAllx library enriched using SureSelect™.  
Analyzed using DNA High Sensitivity Kit.




# PCR amplification – Bead carryover



**Cause** : SPRI bead carryover from post-PCR clean-up step. Use a strong magnet for bead separation and pipette carefully during elution to avoid disturbing beads.


# PCR amplification – Residual Buffer



**Cause** : Likely due to buffer carryover from post-PCR clean-up using columns.




# PCR amplification – Low PCR yield



**Cause** : Inefficient PCR cycling results in low yield and an excess of primers. This can be a result of poor adapter ligation, low DNA quality, inefficient bead clean-up, use of too few cycles or PCR instrument is not well-calibrated. When performing target capture, low yields after post-capture PCR can also indicate suboptimal hybridization.

Based on previous Bioanalyzer runs, it is possible to identify the cause. Ensure that excessive evaporation has not occurred during hybridization (at least 21 $\mu$ l left). If there is sufficient DNA, proceed with sequencing. If not, repeat PCR with remaining library.

# PCR amplification – Chip overloading and spurious peaks



**Cause** : Too much DNA loaded on chip for analysis.

DNA 1000/7500	0.1-50ng/ $\mu$ l
High Sensitivity DNA	100-10000pg/ $\mu$ l

# Amplification-free Library Preparation

Amplification-free Illumina sequencing-library preparation facilitates improved mapping and assembly of (G+C)-biased genomes

Iwanka Kozarewa<sup>1,2</sup>, Zemin Ning<sup>1,2</sup>, Michael A. Quail<sup>1</sup>, Mandy J. Sanders<sup>1</sup>, Matthew Berriman<sup>1</sup> & Daniel J. Turner<sup>1</sup>

- Publication from the Sanger Institute on sequencing of GC-biased genomes using an amplification-free protocol for library construction.
- *P.falciparum* (malaria) and *B.pertussis* genomes are 85% AT-rich and 68% GC-rich, respectively.
- Any PCR-based amplification steps will create significant bias.
- Special adapters were designed to contain complementary sequences to oligos on surface of flow cell.
- Due to the absence of amplification, they used the High Sensitivity DNA Kit to quality control their libraries and QPCR to quantitate.

Library <sup>a</sup>	Organism	Genome size (Mb)	Insert size (bp)	Read length (bp)	Number of reads	Fold raw read coverage	Number of assembled bases
NP-3D7-S	<i>P. falciparum</i> 3D7	23	200	36	28,009,122	43	19,025,823
NP-3D7-L	<i>P. falciparum</i> 3D7	23	200	76	19,556,224	64	21,092,855
STD-PF88	<i>P. falciparum</i> 3D7	23	200	37	110,939,984	174	NA <sup>b</sup>
STD-PF3	<i>P. falciparum</i> 3D7	23	200	37	75,083,768	114	NA <sup>b</sup>
STD-PF2	<i>P. falciparum</i> 3D7	23	200	37	62,802,164	96	NA <sup>b</sup>
STD-PF85	<i>P. falciparum</i> 3D7	23	200	37	13,530,194	21	NA <sup>b</sup>

Kozarewa et al, Nat Methods. 2009 Apr;6(4):291-5.

# Agenda

- 1** Next-gen sequencing library preparation workflow and important quality control steps.
- 2** Introduction to the 2100 Bioanalyzer.
- 3** The use of Bioanalyzer assays for:
  - assessing quality of starting material
  - monitoring size distribution after fragmentation and adapter ligation
  - quantifying yield and detection of artifacts post-PCR amplification
  - detecting small quantities of DNA in amplification-free protocols
- 4** Best practice tips to always obtaining accurate and reproducible data when using the Bioanalyzer.

# Bioanalyzer RNA/DNA Assay Workflow

**Prepare gel-dye mix**




**Prime chip and load samples**



**Load chip and run**



**Quality, Sizing and Quantification data**




# Critical steps to success - Reagents




- All reagents need to be warmed up to room temperature (~ 25°C) for at least 30min prior to preparation. Using cold reagents will affect the run.
- Vortex and spin down gel matrix, dye, marker and ladder to mix well.
- For the RNA assay, both the ladder and RNA samples have to be denatured at 70°C for 2-5min before loading on chip.
- After centrifuging the gel/dye mix, avoid pipetting from the bottom of tube to prevent picking up gel particles. These particles can cause spikes.

**DO NOT use expired reagents!**

# Effects of suboptimal reagent preparation



**'Ghost peaks' due to use of cold reagents.**



**Spikes due to presence of gel particles or use of expired reagents**


# Critical steps to success - Priming



- **Ensure that the syringe clip is in the correct position for your assay eg. highest position for RNA Nano assay and lowest for DNA 1000 assay. Refer to Quick Guides.**
- **Retract plunger to 1ml mark before locking the priming station. Press the silver metal down till you hear a 'Click'. If the plunger is not at 1ml mark, unlock priming station before retracting to prevent drawing up gel/dye.**
- **Upon releasing the syringe clip, the plunger will retract quickly and come to a stop. Count 5sec before pulling plunger back up to 1ml mark.**
- **Do not over- or underprime. Load samples and run chip within 5-10min of priming.**




## Critical steps to success – Maintaining your priming station




- **Replace gasket and clean syringe adapter every 3 months or as needed. The pores should not be clogged with dried gel.**
- **Replace syringe with every new kit.**
- **Perform a pressure test to check condition of priming station. With an unused chip, prime, count 5 sec and release. The plunger should retract to the 0.3ml mark in 1-2 sec.**

# Effects of mispriming



**Wavy and elevated baseline**




**Late migration**

## Critical steps to success – Electrode Cartridge





- The electrode cartridge needs to be completely dry and clean before a run. Liquid or salts on the pins will cause leak currents.
- Prior to running the RNA Nano assay, fill electrode cleaner with 350 $\mu$ l of RNaseZap and soak pins for 30 sec. Clean three times each with 350 $\mu$ l fresh sterile water for 1min and let pins dry. Wash with sterile water once after run. Do not use RNaseZap with the RNA Pico assay unless there is suspected RNase contamination.
- For DNA assays, wash once before and after run with 350 $\mu$ l of sterile water.
- Remove chip within 30 min of run. DO NOT leave chips in Bioanalyzer overnight.
- Do not vortex chips for more than 1min. Check for spills on surface of chip.

# Critical steps to success – Maintaining the electrode cartridge





- Remove pins and clean with sterile water and a soft toothbrush every 3 months or as needed.
- Let pins dry completely before re-assembling the cartridge.
- To check dryness and cleanliness of pins, run a Short Circuit Test using an unused chip. Go to 'Instrument Context', select 'Diagnostics' tab and check 'Short Circuit Test'.

# Effects of a wet or dirty pinset




**Wavy baselines and dips caused by leak currents.**

# Critical steps to success – Data Analysis




- Check that all ladder peaks are present.
- Ensure that Lower and Upper Markers have been assigned correctly.
- Ensure that baselines are flat

# Critical steps to success – Data Analysis



## Critical steps to success – Manually assigning markers



**The Lower Marker is now correctly assigned. Sizing, fragment recognition and quantification are now accurate.**



# Summary

- **The importance of starting a next-gen sequencing experiment with a high-quality library of DNA fragments cannot be over-emphasized.**
- **All major next-gen sequencing and SureSelect™ protocols incorporate multiple quality control steps during library preparation, utilizing Bioanalyzer assays.**
- **RNA and DNA assays are used to :**
  - **Assess quality of starting material**
  - **Monitor size distribution after fragmentation**
  - **Determine success of adapter-ligation**
  - **Determine yield and quality of PCR-amplified libraries**
  - **Detect very small amounts of RNA/DNA in minimally-amplified or amplification-free protocols**
- **It is necessary to adhere to Bioanalyzer assay protocols, as well as to maintain the priming station and electrode cartridge regularly to ensure successful sample analysis and obtain accurate data.**