









Mike Kerber Bio-applications Manager June 2015



## Nikon - Stochastic Optical Reconstruction Microscopy





5nm McGuffee et al. PLoS Comput Biol. 2010.



## 10X Resolution Increase (X, Y, Z) Modalities: 2D, TIRF, 3D, Multi-Color

## How does **STORM** work?



*Emitter isolation:* Only some lights are on at a time
 *Localization:* Pinpoint centroids of visible lights
 *Repeat!*

Requires *photoswitchable* fluorescent probes!

## *"Emitter Isolation"* = Only a subpopulation of fluorophores is on at once

Conventional Fluorescence





Can't see individuals in the crowd...

#### *Emitter Isolation*





#### Individuals distinguishable!

## **Continuous** STORM



Common small molecule dyes will *"blink"* in the presence of *"imaging buffers"* and *high laser powers*.

## With some assumptions, we can accurately find the center of a diffuse spot

 $\otimes$ 

 $\otimes$ 





## Single-Molecule Localization



#### Image of one fluorescent molecule

#### Fitting to a gaussian model



(c)

#### FIONA

Fluorescence Imaging with One Nanometer Accuracy

Paul Selvin



## STORM = *Emitter Isolation* and *Localization*



Rust, Bates & Zhuang, Nat. Methods, 2006 Bates, Huang, Dempsey & Zhuang, Science, 2007



## **Dual Color STORM**

African Green Monkey Kidney (CV-1) Cells with Alexa Fluor 647 (Tubulin) and ATTO 488 (TOMM20)



## **Dual Color STORM**

Human Cervical Cancer Cells (HeLa S3) with Alexa Fluor 647 (NUP153) and ATTO 488 (TPR)





Huang, Zhuang et al, Science (2008)

### 3D Imaging of the Microtubule Network



Scale bar: 200 nm

Huang, Wang, Bates and Zhuang, Science, 2008

5 µm

### Bad Sample = Bad STORM Image

Widefield Clathrin

**Clathrin-STORM** 



Good





## **Determinants of Sample Quality**

1. Probe Choice Dyes that work for N-STORM

2. Labeling Strategies Fixation Immunostaining



\* Don't forget to use glass-bottom (#1.5) dishes to hold STORM imaging buffer!

## Few Dyes Work for N-STORM (Despite what is published)

	Excitation maximum (nm)ª	Emission maximum (nm)ª	Extinction (M <sup>-1</sup> cm <sup>-1</sup> ) <sup>b</sup>	Quantum yield <sup>c</sup>	Detected photons per switching event		Equilibrium on-off duty cycle (400–600 s)		Survival fraction after illumination for 400 s		Number of switching cycles (mean)	
Dye					MEA	βΜΕ	MEA	βΜΕ	MEA	βME	MEA	βME
Blue-absorbing												
Atto 488	501	523	90,000	0.8	1,341	1,110	0.00065	0.0022	0.98	0.99	11	49
Alexa Fluor 488	495	519	71,000	0.92	1,193	427	0.00055	0.0017	0.94	1	16	139
Atto 520	516	538	110,000	0.9	1,231	868	0.0015	0.00061	0.92	0.86	9	17
Fluorescein	494	518	70,000	0.79	1,493	776	0.00032	0.00034	0.51	0.83	4	15
FITC	494	518	70,000	0.8	639	1,086	0.00041	0.00031	0.75	0.9	17	16
Cy2	489	506	150,000	0.12	6,241	4,583	0.00012	0.00045	0.12	0.19	0.4	0.7
Yellow-absorbing												
Су3В	559	570	130,000	0.67	1,365	2,057	0.0003	0.0004	1	0.89	8	5
Alexa Fluor 568	578	603	91,300	0.69	2,826	1,686	0.00058	0.0027	0.58	0.99	7	52
TAMRA	546	575	90,430	0.2	4,884	2,025	0.0017	0.0049	0.85	0.99	10	59
Cy3	550	570	150,000	0.15	11,022	8,158	0.0001	0.0003	0.17	0.55	0.5	1.6
Cy3.5	581	596	150,000	0.15	4,968	8,028	0.0017	0.0005	0.89	0.61	5.7	3.3
Atto 565	563	592	120,000	0.9	19,714	13,294	0.00058	0.00037	0.17	0.26	4	5
Red-absorbing												
Alexa Fluor 647	650	665	239,000	0.33	3,823	5,202	0.0005	0.0012	0.83	0.73	14	26
Cy5	649	670	250,000	0.28	4,254	5,873	0.0004	0.0007	0.75	0.83	10	17
Atto 647	645	669	120,000	0.2	1,526	944	0.0021	0.0016	0.46	0.84	10	24
Atto 647N	644	669	150,000	0.65	3,254	4,433	0.0012	0.0035	0.24	0.65	9	39
Dyomics 654	654	675	220,000	-	3,653	3,014	0.0011	0.0018	0.79	0.64	20	19
Atto 655	663	684	125,000	0.3	1,105	657	0.0006	0.0011	0.65	0.78	17	22
Atto 680	680	700	125,000	0.3	1,656	987	0.0019	0.0024	0.65	0.91	8	27
Cy5.5	675	694	250,000	0.28	5,831	6,337	0.0069	0.0073	0.87	0.85	16	25
NIR-absorbing												
DyLight 750	752	778	220,000	-	712	749	0.0006	0.0002	0.55	0.58	5	6
Cy7	747	776	200,000	0.28	852	997	0.0003	0.0004	0.48	0.49	5	2.6
Alexa Fluor 750	749	775	240,000	0.12	437	703	0.00006	0.0001	0.36	0.68	1.5	6
Atto 740	740	764	120,000	0.1	779	463	0.00047	0.0014	0.31	0.96	3	14
Alexa Fluor 790	785	810	260,000	-	591	740	0.00049	0.0014	0.54	0.62	5	2.7
IRDye 800 CW	778	794	240,000	-	2,753	2,540	0.0018	0.038	0.6	1	3	127

Dempsey et al., 2011

## Where to get secondary antibodies for STORM

cSTORM 2º Antibodies	Sources				
Alexa647	Life Technologies, Jackson				
Cy5	Jackson				
Alexa568	Life Technologies				
Cy3B	DIY (Dye from GE)				
Atto488	Rockland, Sigma				

## Linkage Error: Antibodies separate the fluorophore from the structure of interest



## Our sample prep protocol: a good starting point



## The following *optimizations* are also necessary...

## **Tips** for STORM Sample Preparation

- Compare antibodies from multiple sources.
- *Optimize fixation* (fixative concentration, permeabilization, etc.) to maximize structural preservation and antibody binding.
- *Minimize background* signal levels by titrating primary antibody.
- *Block* with heat-treated sterile filtered blocking serum.
- Don't skip on the *washing steps* and use 1% blocking serum to remove antibodies AT EVERY STEP.
- Lock secondary antibodies in place with *post-staining fixation*.
- Remove residues with *Tween 80 wash*.

## **STORM Sample Prep Summary**

- Use glass-bottom (#1.5) dishes
- Label with Alexa 647 and Atto 488
- Block and wash samples thoroughly
- Use MEA-containing imaging buffer and keep it fresh
- Optimization is usually necessary for the best results!

# TIRF



## Widefield is powerful, but has limitations

#### Out of focus light = high background

#### Low Z resolution

Light is harmful to live cells

## **TIRF** = Total Internal Reflection Fluorescence

#### No out of focus fluorescence

#### Single-molecule sensitivity



#### Plus:

- 150nm Z-resolution
- Reduced phototoxicity to live cells
- Very high speeds possible

So...What's the catch?

### Can only image the 150nm next to the coverslip!

## What near-coverslip applications require *high sensitivity, high speed, and low background?*



#### Focal adhesions



Single molecules



Membrane transport

## How does TIRF work? High angle light reflects off of interfaces

#### Low refractive index

#### **Totally Internally Reflected!**

High refractive index

## TIRF creates an "Evanescent Field" that excites fluorophores within $\sim 150$ nm of the coverslip



#### Coverslip

Oil

**Conventional** Widefield: All fluorophores in the volume are excited Higher signal:noise!





- 1. Focus on your sample interface
- 2. Aim beam directly through sample and onto ceiling
- 3. Focus beam

- This is why you must focus the objective first!
- 4. Adjust beam angle until TIR is achieved



## How do I know if I'm "in TIRF?"

#### Widefield

TIRF



Aperture Plane





What if you want to actually *measure* the penetration depth of the evanescent field?



 $R^2 = X^2 + (R - Z)^2$ 

## Why "High NA?"



🔺 ө → 🦊 Penetration Depth

"Critical angle" = Angle necessary for reflection = sin<sup>-1</sup>(n2/n1) = 63°

> Maximum angle for an objective NA = n sin **e**

"Higher" NA Objective

NA	θ
1.4	67°
1.45	72°
1.49	79°

So, higher NA lenses allow for thinner evanescent fields!

### Also remember: Laser beams have a width



Some laser light at critical angle and some not Total internal reflection!

#### 1.49 NA Objective

*This is the other reason to use "Higher NA" objectives* 

Critical angle threshold

## What if I want to do TIRF with multiple colors?



### Multichannel TIRF considerations

*Remember*: shorter wavelengths of light are refracted more at interfaces.



#### Some options:

- Use *chromatically corrected* optics in your TIRF illuminator
- Use a *motorized* TIRF illuminator
- Use a *separate* TIRF illuminator for each wavelength

Different wavelengths will not strike the back aperture at the same *location* or in the same *focus* 

## What do you need to do TIRF?

- High NA objective
  - 60Х 1.49 Аро
  - 100Х 1.49 Аро
  - 100X 1.49 SR Apo
  - 100X 1.45 K Plan Apo



TIRF cubes



- Laser
- Samples
  - Need a glass-water interface!
  - Get as close to coverslip as possible, minimize ECM!

## Extra hardware for single-molecule live-cell TIRF:

- Maintaining focus is a must!
   Perfect Focus System!
- Single molecules are dim
   EMCCDs
- Molecular motors move fast
  Triggering





- Live-cell fluorophores, like GFP, bleach very fast
  - Minimize exposure to light when not imaging
  - Find focus using some other imaging modality (transmitted, SRIC)

## Common Pitfalls of TIRF

- Very small imperfections in your optics will cause an *interference fringe*
  - DIC prism, etc. can worsen interference patterns and cause ghosting
- Light scatter causes *directional smearing* 
  - Scattering exacerbated by extra coatings
  - Thinner evanescent fields cause less scatter
  - Take images at multiple orientations
- Different wavelengths will not strike the specimen at the same incident angle or same focus
  - Ti TIRF arms have chromatic correction for focus
  - Use different angles for each channel (Motorized or multiple arms)







## **TIRF Summary**



TIRF advantages: high speed, high signal:noise, thin z-sectioning, and low phototoxicity

- *TIRF disadvantages*: can only see ~150nm adjacent to coverslip
- An evanescent field is created by reflecting laser light off of a *glass/water interface*
- Penetration depth depends on wavelength, angle, and refractive index
- Penetration depth can be calculated or measured experimentally
- Objective lenses with at least 1.45 NA are recommended

## Thank you!

#### **Acknowledgements**

Xiaowei Zhuang, Harvard University
 Mike Davidson, FSU
 microscopyu.com

Your local Nikon representatives: ≻Edwin DeFeijter ≻Jim Winkelmann