

Laboratory Diagnostics Core
NE6-217
X51900, x45447
MOUSE CBC/DIFF PROTOCOL

Sample collection

1. Prior to sample collection contact the lab at x51900 to confirm availability of the Advia (our hematology analyzer)
2. A 0.6mL snap cap aliquot tube is ideal for our purposes. It must be a small volume conical tube in order to utilize the entire sample.
3. Add 10 μ L of 0.1M EDTA (18.6g EDTA to 500mL dd water) for each 100 μ L of whole blood.
4. If you are collecting blood from a heart puncture, coat the syringe that you are collecting the whole blood in with the EDTA.
5. If you are collecting with a capillary be sure it has an EDTA coating.
6. Add directly to the aliquot tube that has the EDTA in it. You can start with a smaller amount of EDTA and when you know the blood volume of your collection you can add more to make it 10 μ L EDTA per every 100 μ L whole blood. Immediately **MIX WELL***
7. If you don't have enough volume and you need to dilute your whole blood to get enough volume, dilute in 3% BSA (bovine serum albumin) in PBS prior to analyzing. Mix well.
 - a. Examples:
 - i. x4 dilution 88 μ L of EDTA blood to 264 μ L diluent
 - ii. x8 dilution 44 μ L of EDTA blood to 308 μ L diluent
 - iii. x10 dilution 35 μ L of EDTA blood to 315 μ L diluent
8. The Advia needs 200 μ L for one analysis, 400 μ L of sample will allow for a repeat if needed.
9. Samples should be run within one hour of collection or refrigerated for up to 8 hours.
10. The Advia does not account for your dilution. You will make your calculations.
11. Hours of operation: 7:30am-3:30pm, Monday thru Friday.
12. \$15.00 per sample.
13. Please bring the activity number that you would like billed.

*The most common problem is sample forming small clots before mixing. Your sample can not be evaluated if this happens. You must mix well immediately.