

# Fluorescent Protein (GFP, mCherry etc) Cell Sorting

## Materials

Negative control cells with no fluorescent protein

Presorted cells filtered through 40  $\mu\text{M}$  mesh and resuspended at  $1 \times 10^7$  cells/mL in D-PBS or Sort Staining Buffer in 12x75mm polystyrene tube

### **Sort Staining Buffer**

500ml D-PBS or HBSS ( $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  free)

100U/ml DNase (RNase, protease free) 1 vial used/500ml

7.5ml HEPES (15mM)

16.5ml 30% solution BSA (final concentration 1%)

Penn/Strep (5ml)

2mM EDTA (2ml)

0.2 $\mu\text{m}$  filtered

## Procedure

1. Wash cells out of growth media and resuspend in D-PBS or Sort Staining Buffer at  $2 \times 10^7$  cells/mL
2. Filter cells through 40  $\mu\text{M}$  nylon mesh (can be obtained in the FACS Core)
3. Bring post sort tubes (12 x 75 mm) with the media or buffer you prefer.
  - a. For cell growth following the sort, FBS is recommended
  - b. For RNA or DNA extraction following the sort, D-PBS is recommended