

Detection of the 5' insertion of the ROSA targeting cassette into the ROSA26 locus.

The group (V.T. Chu and R. Kuhn) published a paper entitled “*Efficient generation of Rosa26 knock-in mice using CRISPR/Cas9 in C57Bl/6 zygotes*” in BMC Biotechnology (2016). The primer pairs used to detect insertion of the Rosa cassette into the Rosa26 locus are described.

1. Rosa HR detection:

The primers are:

Rosa26 Fw Out 1.1HA	5' -TAGGGGATCGGGACTCTGGC-3'
Rosa26 Rev CAG	5'- GGGCTATGAACTAATGACCCCG-3'
Rosa26-Fw bGH	5'- GGAAATTGCATCGCATTGTCTGA-3'
Rosa26-Rev Out 4.3HA	5'- ACCATTCTCAGTGGCTCAACAACA-3'

PCR conditions:

Input-

Genomic DNA	2.0ul
2x GC buffer I	25ul
dNTPs (2.5mM)	8ul
Primer F (100uM)	.25ul
Primer R (100uM)	.25ul
LA Taq	0.5ul
H ₂ O	14.0ul
total volume	<u>50.0ul</u>

Program:

94°C	1 min	
94°C	30 sec	} 30 cycles
60°C	30 sec	
72°C	2 min	

72°C 7 min

4 Hold

2. Amp^r random insertion detection:

The integration of the Rosa targeting cassette may have occurred randomly, into the desired location in the *Rosa26* locus or you may have both types of integration events in the same genome. To identify random integrations, you can screen for the presence of the Amp^r cassette as the vector DNA would have integrated along with the targeting cassette. Use the same **LA Taq DNA polymerase** (Takara # RR02AG) for the analysis and use Amp^r specific primers.

Amp^r Primers:

Amp ^r (A) Forward	5'-AAGATGCTGAAGATCAGTTG-3'	>143bp
Amp ^r (B) Reverse	5'-ATAATACCGCGCCACATAGC-3'	
Amp ^r (C) Forward	5'-AATTAATAGACTGGATGGAG-3'	>199bp
Amp ^r (D) Reverse	5'-TTCATCCATAGTTGCCTGACTC-3'	

Program:

95°C	10 min	} 22 cycles
94°C	1 min	
53°C	1 min	
72°C	1 min	
72°C	10 min	
4°C	HOLD	

PCR #_390__

PCR DETECTION OF TRANSGENES (Genotyping)

Date_10/11/2022__

I. Insertion: Rosa26-Color Switch -- To test specific insertion in Rs26

II. Primers: *Fw1 OUT Rs26* (5'-taggggatcgggactctggc-3')
Rev1 CAG (5'-gggctatgaactaatgaccccg..3')

IV. MIX: LA Taq DNA w/ GC buffer master mix 48uL each + 2uL Sample

<u>x 1 Sample</u>		<u>x5 Samples</u>
<i>Primer 1 (100uM)</i>	0.25uL	1.25
<i>Primer 2 (100uM)</i>	0.25uL	1.25
LA Taq	0.5uL	2.5
2x GC buffer I or II	25uL	125
dNTP (2.5mM)	8uL	40uL
H2O	15uL	75uL

V. PCR program:

94°C 1min
94°C 30sec
60°C 30sec
72°C 2min
72° 5min
4 Hold

30cycles


VI. Samples:

1. __927__
2. __934__
3. __936__
4. __H2O__
5. __939 (WT)__
6. _____

III. Results:

WT band: N/A

KI band (R26 Specific): 1248bp



Buffer I: Fragments longer than 5Kb
Buffer II: Fragments 2-3kb in length

PCR #_390__

PCR DETECTION OF TRANSGENES (Genotyping)

Date_10/11/2022__

I. Insertion: Rosa26-Color Switch -- To test specific insertion in Rs26

II. Primers: *Fw1 OUT Rs26* (5'-taggggatcgggactctggc-3')
Rev1 CAG (5'-gggctatgaactaatgaccccg..3')

III. Results:

WT band: N/A

KI band (R26 Specific): 1248bp

IV. MIX: KOD hot start master mix 18uLeach + 2uL Sample

<u>x 1 Sample</u>		<u>x5 Samples</u>
(100uM) <i>Primer 1</i>	0.25uL	<u>1.25</u>
(100uM) <i>Primer 2</i>	0.25uL	<u>1.25</u>
<i>MM</i>	10uL	<u>50</u>
<i>H2O</i>	8.5uL	<u>42.5</u>

V. PCR program:

95°C 2min
95°C 20sec
60°C 10sec
70°C 1.15min
70° 2min
4 Hold

30cycles

VI. Samples:

1. __927__
2. __934__
3. __936__
4. __H2O__
5. __939 (WT)__
6. _____

