

Detection of the 5' and 3' insertion of the ROSA targeting cassette into the ROSA26 locus.

The group (V.T. Chu and R. Kuehn) published a paper entitled "*Efficient generation of Rosa26 knock-in mice using CRISPR/Cas9 in C57Bl/6 zygotes*" in BMC Biotechnology (2016). The primer pairs used to detect insertion of the Rosa cassette into the Rosa26 locus are described.

***Please note the difference in the 5' and 3' PCR Taq kits

For 5' Detection use Takara LA Taq with GC Buffer (#RR02AG)

The 5' primers are:

Rosa26 Fw Out 1.1HA	5' -TAGGGGATCGGGACTCTGGC-3'
Rosa26 Rev CMV	5'- GGGCTATGAACTAATGACCCCG-3'

5' PCR conditions:

Input-	
Genomic DNA	2.0ul
2x GC buffer I or II	25ul
dNTPs (2.5mM)	8ul
Primer F (100uM)	.25ul (or 5ul @10pmol/ul)
Primer R (100uM)	.25ul (or 5ul @10pmol/ul)
LA Taq	0.5ul
H ₂ O	14.0ul
total volume	<u>50.0ul</u>

Program:

94°C	1 min	
98°C	30 sec	} 30 cycles
68°C	30 sec	
72°C	2 min	

72°C 7 min
4 Hold

For 3' Detection use Takara LA Taq DNA Polymerase (#RR002A)

The 3' primers are:

Rosa26-Fw bGH	5'- GGAAATTGCATCGCATTGTCTGA-3'
Rosa26-Rev Out 4.3HA	5'- ACCATTCTCAGTGGCTCAACAACA-3'

3' PCR conditions:

Input-

Genomic DNA	2.0ul	
10x LA Buffer II	5ul	
25mM MgCl ₂	5ul	
dNTPs (2.5mM)	8ul	
Primer F (100uM)	.25ul	(or 5ul @10pmol/ul)
Primer R (100uM)	.25ul	(or 5ul @10pmol/ul)
LA Taq	0.5ul	
H ₂ O	29.0ul	
total volume	50.0ul	

Program:

94°C 1 min

98°C	10 sec	] 30 cycles
68°C	15 min	

72°C 7 min

4 Hold

The integration of the Rosa targeting cassette may have occurred randomly, into the desired location in the *Rosa26* locus or you may have both types of integration events in the same genome. To identify random integrations, you can screen for the presence of the Amp cassette as the vector DNA would have integrated along with the targeting cassette. Use the same **LA Taq with GC Buffer** (Takara # RR02AG) for the analysis and use AMP specific primers.

Amp^r Primers:

Amp ^r (A) Forward	5'-AAGATGCTGAAGATCAGTTG-3'	>143bp
Amp ^r (B) Reverse	5'-ATAATACCGCGCCACATAGC-3'	
Amp ^r (C) Forward	5'-AATTAATAGACTGGATGGAG-3'	>199bp
Amp ^r (D) Reverse	5'-TTCATCCATAGTTGCCTGACTC-3'	

Program:

95°C	10 min	} 22
94°C	1 min	
53°C	1 min	
72°C	1 min	
72°C	10 min	
4°C	HOLD	

PCR DETECTION OF TRANSGENES II

Date: 11/25/25

I. Transgene: Test for 5' Rosa HR insertion
 Primers: R26 FW Out 1.1HA, R26 REV-CMV

V. Results: photo 1.4% gel

II. PCR Rxn. (Takara LA Taq w/ GC Buffer) 8 X


DNA	_1_ μ l		
2x GC Buffer II	12.5 μ l		100
dNTPs	4.0 μ l		32
Olgo A <u>F</u>	2.5 μ l		20
Olgo B <u>R</u>	2.5 μ l		20
Olgo C _____	_____ μ l		_____
Olgo D _____	_____ μ l		_____
Ster. H ₂ O	2.25 μ l		18
Taq. Pol	0.25 μ l		2

(25 μ l / rxn)

III. PCR program: Standard Biometra #3,5
 Other: 60°C (30s, 30s, 2min) _____
30 cycles _____

IV. Samples:

Lane	1. <u>H3/RI</u>	8. <u>1Kb Direct Load</u>	15. _____	22. _____	29. _____	36. _____
	2. <u>54</u>	9. _____	16. _____	23. _____	30. _____	37. _____
	3. <u>64</u>	10. _____	17. _____	24. _____	31. _____	38. _____
	4. <u>65</u>	11. _____	18. _____	25. _____	32. _____	39. _____
	5. <u>76</u>	12. _____	19. _____	26. _____	33. _____	40. _____
	6. <u>NC</u>	13. _____	20. _____	27. _____	34. _____	
	7. <u>H2O</u>	14. _____	21. _____	28. _____	35. _____	



PCR DETECTION OF TRANSGENES II

Date: 11/19/2025

I. Transgene: Test for 3' Rosa HR insertion


Primers: R26 bGH-F, R26 Out 4.3HA-R

V. Results: photo

II. PCR Rxn. (Takara LA Taq DNA Poly) 8 x

DNA	<u>1</u> μ l	<u> </u>
10x LA Buffer II	<u>2.5</u> μ l	<u>20</u>
25mM MgCl ₂	<u>5</u> μ l	<u>20</u>
dNTPs	<u>0.5</u> μ l	<u>32</u>
Oligo A <u> F </u>	<u>2.5</u> μ l	<u>20</u>
Oligo B <u> R </u>	<u>2.5</u> μ l	<u>20</u>
Oligo C <u> </u>	<u> </u> μ l	<u> </u>
Oligo D <u> </u>	<u> </u> μ l	<u> </u>
Ster. H ₂ O	<u>9.75</u> μ l	<u>78</u>
Taq. Pol	<u>0.25</u> μ l	<u>2</u>

(25 μ l / rxn)



III. PCR program: Standard Biometra #3,10

Other: 68°C (10sec, 15min)

30 cycles

3'Rosa HR 4662bp

IV. Samples:

Lane1.	<u>H3/RI</u>	8.	<u>1Kb Plus</u>	15.	<u> </u>	22.	<u> </u>	29.	<u> </u>	36.	<u> </u>
2.	<u>54</u>	9.	<u> </u>	16.	<u> </u>	23.	<u> </u>	30.	<u> </u>	37.	<u> </u>
3.	<u>64</u>	10.	<u> </u>	17.	<u> </u>	24.	<u> </u>	31.	<u> </u>	38.	<u> </u>
4.	<u>65</u>	11.	<u> </u>	18.	<u> </u>	25.	<u> </u>	32.	<u> </u>	39.	<u> </u>
5.	<u>76</u>	12.	<u> </u>	19.	<u> </u>	26.	<u> </u>	33.	<u> </u>	40.	<u> </u>
6.	<u>NC</u>	13.	<u> </u>	20.	<u> </u>	27.	<u> </u>	34.	<u> </u>		
7.	<u>H2O</u>	14.	<u> </u>	21.	<u> </u>	28.	<u> </u>	35.	<u> </u>		

PCR DETECTION OF TRANSGENES II

Date: 1/6/25

I. Transgene: Amp Detection for Rosa 26
 Primers: Amp^r -A (Forward), B (Reverse)

V. Results: photo 2% gel

II. PCR Rxn. (Takara LA Taq w/ GC Buffer) 4 X

DNA	<u>1</u> μ l	<u> </u>
2x GCBuffer II	<u>12.5</u> μ l	<u>50</u>
dNTPs	<u>4.0</u> μ l	<u>16</u>
Oligo A <u>A</u>	<u>2.5</u> μ l	<u>10</u>
Oligo B <u>B</u>	<u>2.5</u> μ l	<u>10</u>
Oligo C <u> </u>	<u> </u> μ l	<u> </u>
Oligo D <u> </u>	<u> </u> μ l	<u> </u>
Ster. H ₂ O	<u>2.25</u> μ l	<u>9</u>
Taq. Pol	<u>0.25</u> μ l	<u>1</u>

(25 μ l / rxn)

III. PCR program: Standard Biometra #64

Other: 53°C (1,1,1)

22 cycles

143bp

IV. Samples:

Lane	1. <u>100bp</u>	8. <u> </u>	15. <u> </u>	22. <u> </u>	29. <u> </u>	36. <u> </u>
	2. <u>NC (C57Bl6)</u>	9. <u> </u>	16. <u> </u>	23. <u> </u>	30. <u> </u>	37. <u> </u>
	3. <u>PC (1ng/μl pR26)</u>	10. <u> </u>	17. <u> </u>	24. <u> </u>	31. <u> </u>	38. <u> </u>
	4. <u>H2O</u>	11. <u> </u>	18. <u> </u>	25. <u> </u>	32. <u> </u>	39. <u> </u>
	5. <u> </u>	12. <u> </u>	19. <u> </u>	26. <u> </u>	33. <u> </u>	40. <u> </u>
	6. <u> </u>	13. <u> </u>	20. <u> </u>	27. <u> </u>	34. <u> </u>	
	7. <u> </u>	14. <u> </u>	21. <u> </u>	28. <u> </u>	35. <u> </u>	

