

Running a microCT Scan (for zebrafish)

Part 1: microCT Preparation

1. Turn on the microCT computer (the black console, then the monitor) and start the GPUReconServer background server. Then start the SkyScan 1275 program. Check that the microCT hard drive is plugged into the USB port on the left side of the monitor.
2. Click the X-Ray icon ☢️ to warm up the scanner. This could take anywhere from 15 minutes to over an hour, depending on when the scanner was last used, so plan accordingly.

Part 2: Sample Preparation

1. Refer to [PTA Staining Protocol](#). Samples should be fixed for at least 24 hours in ethanol in the 4°C fridge prior to scanning.
 - a. To fix fish, euthanize them with Tricaine before transferring them to PFA for 1-2 days.
 - b. Then, replace the PFA with increasing graded concentrations of ethanol, 30%-50%-70% for 2 hours each, before placing the fish in 0.7% PTA in 90% methanol for a week.
 - c. Wash once with ethanol and discard. The fish can be stored in ethanol indefinitely.
2. Select a container that best fits the size of the sample:
 - a. 1.5 mL Eppendorf tube: 8.0 um resolution
 - b. 15 mL conical tube: 10.5 um resolution
 - c. 50 mL conical tube: 35.5 um resolution
 - d. Gelato jar: 48.5 um resolution
3. Dry the sample **thoroughly** with Kimwipes and place the sample in purple low-density foam. The sample can also be wrapped in Kimwipes. Note the position of the head and tail.
 - a. If wrapping the sample in Kimwipes:
 6. Dry the fish superficially on paper. Then wrap the fish in a roll of tissue paper, with a width cut to about 50% longer than the fish, as shown in figure 3.



Figure 3. After removal from storage, the fish is laid on a strip of tissue paper and rolled up prior to placing in a tube sample holder for microCT imaging.

i.

4. Place the sample tail-first into the tube so that when the tube is screwed into the base, the head is pointing downwards.
5. Wet the cotton gauze at the bottom of the tube with 1-2 drops of PBS.
6. Cover the top of the test tube (if exposed) with parafilm.
7. Mount the sample onto the imaging stage. Adjust the sample so that the tube is as straight as possible. Rotate the stage to check from every angle.

Part 3: Imaging Preparation

1. Turn off the X-Ray  when finished warming up. Turn on the X-Ray  and Grab  functions.
2. Perform a flat-field correction.
3. Type **Ctrl + Shift + Alt + S** to unlock scanning settings. Go to Options → Scanning Settings and adjust as needed, only changing the values in the high resolution column.
 - a. Generally for zebrafish:
 - i. **45 kV voltage**
 - ii. **222 uA current**
 - iii. **0.2° rotation step**
4. Click “Ok” to save the scanning settings.
5. Double-click ff in the upper left corner to turn flat-field correction off. Turn the X-Ray  on and right-click in the middle of the screen. The average exposure should be between 61-62%.
6. Click “Update flat-field references” and select the first and last boxes. When the computer is finished updating, turn flat field correction back on. The average exposure should be between 90-92%.
7. Turn off the X-Ray  and open the chamber  to place the sample back in. Turn on the X-Ray  and Grab  functions.
8. Move the sample using the up and down arrows, and change the size of the sample by adjusting the resolution.
9. Rotate the sample to check for white space (air) on both sides of the tube. The top part of the tube will be more affected by any crookedness or wobbling from the base.
10. Set the upper and lower bounds of the scan by double-clicking a little above and below the fish.
11. Press Start Scan  and a directory for naming the file should show up.
 - a. The name should include:
 - i. A brief description listing any experimental conditions
 - ii. The date of the scan
 - iii. The SL of the fish measured from snout to base of tail (if applicable)
 - b. Example: PracticeScan_4April25_10.5um_45kV_222uA in Tho’s folder
12. Click “Make new folder” with the file name, and then click out of the folder onto a random file before clicking back in. The folder name should update automatically (the computer is old and laggy, so this is a safeguard to make sure the file is named properly).
13. Copy-paste the new name into the directory to double-check the file name.
14. Check the first four boxes in the Start Scan box. Then, hit Scan to start.
 - a. Scans can take anywhere from 45 minutes to 2 hours, depending on the size of the fish and the region to be scanned (whole body, head only, or tail only).
 - b. Plan scans accordingly and set a timer to turn off the scanner and remove the sample from the chamber. Alternatively, click “Turn computer off after finished scanning.”

Part 4: Reconstruction

1. Go to File Explorer and select any one of the microCT images from the desired scan.
2. Turn on the NRecon software. The GPUReconServer should run automatically in the background.
3. Go to the **Start tab** → **Preview fastest slice**. This selects one representative 2D cross-section from the raw scanning data to be processed so you can more quickly see the effects of various transformations.
4. Go to **View** → **Show Profile Window**, which should show a graph with black and green lines.
 - a. If the black and green lines match up, no action needed.
 - b. If the graph is a complete mess, rescan the fish.
5. Perform an X/Y correction. Click **Actions** → **X/Y alignment with reference scan** → adjust the box to fit the volume of the tube. Click “Match”, wait for the system to finish, then click “Accept”.
6. Click **Start** → **Preview** → **Output tab** and right-click to draw a line bisecting the fish. Then click Update Rotation.
 - a. **“X/Y alignment”** (misalignment compensation) is a global, uniform shift of the sample’s rotation axis relative to the detector. This is needed when there are minor inaccuracies in scanner calibration or if the sample was mounted off-center.
 - b. **“Update rotation”** addresses a non-uniform misalignment that varies along the Z-axis of the sample. This correction is needed when the sample’s rotation axis is not perfectly perpendicular to the detector plane (like the Leaning Tower of Pisa), so that features are more blurred at the top and/or bottom of the sample.
 - i. We usually do both corrections for zebrafish samples, but usually only X/Y alignment is needed.
7. Go to **Start** and adjust the green slider to the widest part of the sample.
 - a. For the posterior (tail) end of zebrafish, this is the section closest to the peduncle.
 - b. For the anterior (head) end of zebrafish, aim for the ear bones (otoliths), which should appear as one or two very dense structures made of calcium carbonate.
 - c. For samples that have a relatively constant width and density, move the green slider to a representative cross-section.
8. Click “Fine Tuning”, select “Post-alignment”, and choose 5 trials with a parameter step of 0.5.
9. Go to **Output** and resize the blue circle so that it fits inside the tube. You can increase contrast by moving the right slider on the histogram inwards.
 - a. The x-axis on the histogram represents the X-ray attenuation coefficient, how effectively a material absorbs or scatters x-rays as they pass through it. Denser materials absorb more x-rays and have higher attenuation coefficients, while less dense materials will have lower attenuation coefficients.
 - b. The y-axis represents the number of pixels that have a particular intensity value.
 - c. **Don’t adjust the left slider!**
10. After resizing the blue circle, move the right slider back to roughly its original position (make sure the right slider is completely outside of the black lines).
11. Use the left and right sliders at the top of the screen to click back and forth. Select the clearest image with the highest resolution, then leave the selected image open on the screen.

12. Return to the **Start** tab and repeat the **Fine Tuning** → **Post-alignment** process for each section.
13. When done choosing the clearest slices, adjust the red bounds on the Start tab so that they are closer to the top and bottom of the sample, leaving roughly 1-2 cm of buffer.
14. On the **Start** tab, click “Start” to begin the reconstruction process. The reconstruction data will be stored in the same folder as the scan data, in a new folder called “ScanTitle_Rec”.

Part 5: 3D Visualization and Analysis

Introduction

There are many, many different ways to visualize and analyze the reconstructed results of microCT scans. In the past, the lab has used Amira-Avizo software from Thermo-Fisher to create 3D renderings of microCT scans – however, technical difficulties accessing the cloud network and slow computer processing speeds due to high volumes of data have made using Amira for scans impractical.

The lab owns a [Bruker SkyScan 1275](#) benchtop microCT scanner. The scanner comes with built-in software such as **NRecon** (reconstructs 2D projection images from a microCT scan into 3D volumes), **DataViewer** (allows for slice-by-slice visualization of 3D volumes), **CTan** (CT-analyzer, useful for density and morphological analysis), **CTvol** (CT-volume, provides surface rendering), and **CTvox** (provides volume rendering).

NOTE: All user manuals and Bruker resources can be found in folders on the microCT computer! Phil Salmon’s tutorial videos on YouTube are also incredibly helpful.


▶ Overview of all software in Bruker MicroCT 3D Suite: what does each program do?

Calculating Bone Mineral Density (BMD) on CTan


▶ CT-Analyser ("CTAn") - an overview and introduction.

▶ Bruker microCT tutorial: Calibration and measurement of bone mineral density (BMD and TMD) in...


1. First, open the CTan software. Go to File → Open... → choose the reconstruction file (.bmp) you want to import. Make sure **Open as Dataset** is selected.
2. You should see a **Shadow Projection** tab showing the 2d projections obtained from the microCT scan, a **Dataset** tab displaying each reconstruction file, and a **Cross Section** tab displaying the reconstructed 3d cross sections.
3. Adjust the slider to the left of the **Dataset** tab to view different cross-sections.
4. In the upper toolbar, click Regions of interest preview  to set your desired region of interest. For example, if scanning a density phantom, adjust the slider until you can see all 5 density phantoms in the same cross-section. Then right-click on the sun icon  → Set the Top/Bottom of Selection to set the upper and lower bounds of the ROI.
5. On the **Cross Section** tab, either left-click and drag to make a shape, or choose from a list of shapes to create the ROI.



- When creating the ROI, select a range well within the outer borders of the material to avoid intermediate areas which may complicate calculations.



Round ROI of a density phantom within a selected range.



The ROI is well contained within the borders of the material.

- Set an ROI at the upper and lower bounds. The software will automatically generate "Interpolated ROIs" to dynamically connect between the endpoints. Double-check the size and range of the ROI with the **Dataset** slider.
- Go to **Binary Selection**  → Histogram tab → From dataset  to generate binary selections using the entire dataset. Toggle VOI view  to make sure the ROI is selected (background surrounding the ROI should be green, not black).
- From here, you can obtain the grayscale indexes, Hounsfield units, attenuation coefficient, and bone mineral density for the ROI. To do this, scroll all the way down until you see the **Mean (total)** value. This is the mean value of your metric of interest.