



University of Pittsburgh Cancer Institute

Lentiviral Core Facility

Recommended Protocol

Director: Robert W. Sobol, PhD
Hillman Cancer Center
5117 Centre Ave., Suite 2.6
Pittsburgh, PA 15213
Phone: (412) 623-7764 / Fax: (412) 623-7761
Email: rws9@pitt.edu

Recommended Protocol

HIV Lentivirus Transduction (suspension cells)

Summary: HIV-based lentivirus particles are prepared by transfection of 4 plasmids into 293-FT cells [Dull T et al. "A third-generation lentivirus vector with a conditional packaging system" J Virol. 1998 Nov;72(11):8463-71 and Robinson DA, Dillon CP, Kwiatkowski AV, et al. A lentivirus-based system to functionally silence genes in primary mammalian cells, stem cells and transgenic mice by RNA interference" Nature Genetics 2003;33(3):401-6].

Virus particles are isolated from the cell culture supernatant and may be stored at -80°C or used immediately for transduction of the target cell. As per current University of Pittsburgh IBC/rDNA protocols, this virus preparation may be used following BSL-2 bio-safety procedures.

Day 1:

1) Prepare the following dilutions on ice in sterile tubes:

<u>Well #</u>	<u>Vol. DMEM*</u> (ml)	<u>Vol. Virus</u> (ml)	<u>Polybrene</u> (µl)
1 (negative control)	2.0	0	2
2	1.0	1.0	2

(The concentration of the stock solution of polybrene is 8 mg/ml, which is diluted down to 8 µg/ml when added to the virus. For example, 3 µls of the polybrene is added to 3 mls of the Lentivirus stock)

*or you may use your normal cell growth media.

2) Carefully pellet 125,000 cells in a 15 ml tube (1500xg, 5 min). Remove media and suspend cells with the virus solution (2 ml) as indicated in the table. Transfer to a 60 mm dish.

3) Incubate at 32°C/ 5% CO₂ overnight.

Day 2:

Carefully pellet 125,000 cells in a 15 ml tube (1500xg, 5 min). Remove media and suspend cells with fresh growth media, transfer to a 60 mm dish. **Now incubate at 37°C.**

Days 5-8: Analyze cells for GFP expression on microscope or protein expression by extract preparation and immunoblot.

It may take at least 48-96 hrs for transgene to express.