

BD- C6 analyzer

Filter options for alternate laser configurations

Standard: 3-blue 1-red

Detector Position	Laser ₁	Filter Options ₂	Example Fluorophores
FL1	488 nm	533/30 or 530/30	FITC, GFP, YFP, CFSE, Alexa Fluor® 488
FL2	488 nm	585/40	PE, PI, PE-CF594
FL3	488 nm	670 LP	PerCP, PE-Cy5, PerCP-Cy5.5, PE-Cy7, PI, 7-AAD
FL4	640 nm	675/25	APC, Alexa Fluor® 647