

FACScan Configuration			
Laser	Bandpass filter	Detector Name	*Suggested Fluorochromes
Blue (488nm)	488/10	FSC	
	488/10	SSC	
	530/30	FL 1	Fitc, Alexa488, GFP, YFP
	585/42	FL 2	PE, PI
	LP670	FL 3	PerCP, PerCP Cy5.5, PE Cy5, 7AAD, PE Cy7
Red (637nm)	670/20	FL4	APC, Alexa647, Alexa633, To-Pro-3
*List of suggested Fluorochrome names is not exhaustive, please consult FCCF staff for information on other usable reagents.			