

LSRII

LASER	PMT#	DICHROIC FILTER	EM RANGE
BLUE (Ar488)	A	670LP 695/40	670-715nm
	B	505LP 515/20	505-525nm
	C	- 488/10	483-493nm
RED (639 HeNe)	A	755LP 780/60	755-810nm
	B	690LP 730/45	690-752nm
	C	- 660/20	650-670nm
GREEN (532)	A	735LP 780/60	750-810nm
	B	600LP 610/20	600-620nm
	C	- 575/26	562-588nm
UV (355 HeCd)	A	505LP 525/50	505-550nm
	B	- 450/50	425-475nm
	C	- -	-
Violet (405)	A	535LP 560/40	535-580nm
	B	- 450/50	425-475nm
	C	- -	-

**Please consult the Flow Core or
<https://fluorofinder.com/> to design your
panel using Core specific instruments.**