

Melody

LASER	PMT#	DICHROIC FILTER	EM RANGE
BLUE (Ar488)	A	665LP 700/54	673-727nm
	B	507LP 527/32	511-543nm
	C		
RED (640 HeNe)	A	752LP 783/56	755-811nm
	B	- 660/10	655-665nm
GREEN (532)	A	752LP 783/56	755-811nm
	B	665LP 697/58	668-726nm
	C	605LP 613/18	604-622nm
	D	582/15	575-589nm

**Please consult the Flow Core or
<https://fluorofinder.com/> to design your
panel using Core specific instruments.**