

## Compensation Beads

**Note: It is important to know the species and isotype of your antibody prior to staining with compensation beads.**

### **For Mouse Isotype Antibodies**

COMPtrol Kit (Spherotech Cat# CMlgP-30-2K)  
Goat anti-Mouse Ig (H+L)  
**Kathy's recommended beads**

CompBeads, Mouse Ig (BD Biosciences # 552843)  
Anti-Mouse Ig, kappa light chain

AbC™ Anti-Mouse Bead Kit (ThermoFisher # A10344)  
Anti-Mouse Ig

### **For Multiple Species Isotype Antibodies**

CompBeads Rat and Hamster Ig, (BD Biosciences Cat# 552845)  
Anti-Rat, Hamster, Ig, kappa light chain  
**Kathy's recommended beads**

VersaComp (Beckman Coulter Cat# B22804)  
Anti-Mouse, Rat, Hamster

UltraComp eBeads™ (ThermoFisher Cat# 01-2222-42)  
Anti-Mouse, Rat, Hamster  
**Note: Low bead concentration. Not recommended for inexperienced users**

OneComp eBeads™ (ThermoFisher Cat# 01-1111-42)  
Anti-Mouse, Rat, Hamster

UltraComp eBeads™ Plus (ThermoFisher Cat# 01-3333-42)  
Anti-Mouse, Rat, Hamster, Rabbit, Human  
**Note: Low bead concentration. Not recommended for inexperienced users**

AbC™ Total Antibody Compensation Bead Kit (ThermoFisher # A10513)  
Anti-Mouse, Rat, Hamster, Rabbit

AbC™ Anti-Rat/Hamster Bead Kit (ThermoFisher # A10389)  
Anti-Rat, Hamster

### **Specialty Compensation Beads**

ArC™ Amine Reactive Beads (ThermoFisher Cat# A10628)  
Live/Dead Fixable Viability Dye  
**Note: You will need to stain the ArC beads with more dye than your cells**

GFP BrightComp eBeads (ThermoFisher Cat# A10514)  
GFP expressing cells